



# Validation Data Sheet

dPCR CNV Probe Assay	NRAS
Assay Catalog No.	DCG0000201
dPCR Instrument	QIAcuity digital PCR instrument
dPCR Master Mix	QIAcuity Probe PCR Master Mix
Experimental Sample(s)	genomic DNA, Xbal-digested

**1D Scatter Plot:** Single-well data of target assay NRAS (Green Channel) and reference assay TERT (Yellow Channel) in different gDNA input amounts (A10,B10,C10: 6ng/rxn, D10,E10,F10: 4ng/rxn, G10,H10: NTCs)



**2D Scatter Plot:** Single-well data of target assay NRAS (Green Channel, X-axis) and reference assay TERT (Yellow Channel, Y-axis) in 4ng gDNA input amount.

