


# Validation Data Sheet

<b>dPCR CNV Probe Assay</b>	<b>EML4-AS1</b>
Assay Catalog No.	DCG0000548
dPCR Instrument	QIAcuity digital PCR instrument
dPCR Master Mix	QIAcuity Probe PCR Master Mix
Experimental Sample(s)	genomic DNA, Xbal-digested

**1D Scatter Plot:** Single-well data of target assay EML4-AS1 (Green Channel) and reference assay TERT (Yellow Channel) in different gDNA input amounts (A9,B9,C9: 4ng/rxn, D9,E9,F9: 6ng/rxn, G9,H9: NTCs)



**2D Scatter Plot:** Single-well data of target assay EML4-AS1 (Green Channel, X-axis) and reference assay TERT (Yellow Channel, Y-axis) in 4ng gDNA input amount.

