

# RT<sup>2</sup> Profiler PCR Array (Rotor-Gene<sup>®</sup> Format)

## Rhesus Macaque Housekeeping Genes

Cat. no. 330231 PAQQ-000ZR

For pathway expression analysis

Format	For use with the following real-time cyclers
RT <sup>2</sup> Profiler PCR Array, Format R	Rotor-Gene Q, other Rotor-Gene cyclers

### Description

The Rhesus Macaque Housekeeping Genes RT<sup>2</sup> Profiler PCR Array profiles the expression of 8 identical sets of 12 commonly used housekeeping genes. This array can be easily used to identify genes with a constant level of expression among different experimental conditions for use in normalizing relative gene expression profiling experiment. Housekeeping genes encode proteins that are usually essential for the maintenance of cellular function. Their expression often, but not always, remains constant under most experimental conditions. Some of these genes may be regulated under certain circumstances and may vary with cell type. Using real-time PCR, research studies can easily and reliably analyze the expression of these housekeeping genes among different experimental conditions with this array.

For further details, consult the *RT<sup>2</sup> Profiler PCR Array Handbook*.

### Shipping and storage

RT<sup>2</sup> Profiler PCR Arrays in the Rotor-Gene format are shipped at ambient temperature, on dry ice, or blue ice packs depending on destination and accompanying products.

For long term storage, keep plates at –20°C.

**Note:** Ensure that you have the correct RT<sup>2</sup> Profiler PCR Array format for your real-time cycler (see table above).

**Note:** Open the package and store the products appropriately immediately on receipt.



## Array layout

The 96 real-time assays in the Rotor-Gene format are located in wells 1–96 of the Rotor-Disc™ (plate A1–A12=Rotor-Disc 1–12, plate B1–B12=Rotor-Disc 13–24, etc.). To maintain data analysis compatibility, wells 97–100 do not contain real-time assays but will contain master mix to account for weight balance.

## Gene table: RT<sup>2</sup> Profiler PCR Array

Position	UniGene	GenBank	Symbol	Description
A01	Mmu.4974	NM_001033084	ACTB	Actin, beta
A02	Mmu.5037	NM_001047137	B2M	Beta-2-microglobulin
A03	Mmu.3145	XM_001105471	GAPDH	Glyceraldehyde-3-phosphate dehydrogenase
A04	Mmu.9544	XM_001087699	GUSBB	Glucuronidase, beta, b
A05	Mmu.142	NM_001195533	HSP90AB1	Heat shock 90kDa protein 1, beta
A06	Mmu.12316	XM_001097691	LOC709186	Hypoxanthine-guanine phosphoribosyltransferase-like
A07	Mmu.4112	XM_001102090	LOC712934	Hypothetical LOC712934
A08	Mmu.2957	XM_001091198	NONO	Non-POU domain containing, octamer-binding
A09	Mmu.4126	NM_001195658	PGK1	Phosphoglycerate kinase 1
A10	Mmu.17244	XM_001087648	PPIH	Peptidylprolyl isomerase H (cyclophilin H)
A11	Mmu.2512	XM_001115079	RPL13A	Ribosomal protein L13A
A12	Mmu.12238	XM_001085368	TBP	TATA box binding protein
B01	Mmu.4974	NM_001033084	ACTB	Actin, beta
B02	Mmu.5037	NM_001047137	B2M	Beta-2-microglobulin
B03	Mmu.3145	XM_001105471	GAPDH	Glyceraldehyde-3-phosphate dehydrogenase
B04	Mmu.9544	XM_001087699	GUSBB	Glucuronidase, beta, b
B05	Mmu.142	NM_001195533	HSP90AB1	Heat shock 90kDa protein 1, beta
B06	Mmu.12316	XM_001097691	LOC709186	Hypoxanthine-guanine phosphoribosyltransferase-like
B07	Mmu.4112	XM_001102090	LOC712934	Hypothetical LOC712934
B08	Mmu.2957	XM_001091198	NONO	Non-POU domain containing, octamer-binding
B09	Mmu.4126	NM_001195658	PGK1	Phosphoglycerate kinase 1
B10	Mmu.17244	XM_001087648	PPIH	Peptidylprolyl isomerase H (cyclophilin H)
B11	Mmu.2512	XM_001115079	RPL13A	Ribosomal protein L13A
B12	Mmu.12238	XM_001085368	TBP	TATA box binding protein
C01	Mmu.4974	NM_001033084	ACTB	Actin, beta
C02	Mmu.5037	NM_001047137	B2M	Beta-2-microglobulin
C03	Mmu.3145	XM_001105471	GAPDH	Glyceraldehyde-3-phosphate dehydrogenase
C04	Mmu.9544	XM_001087699	GUSBB	Glucuronidase, beta, b
C05	Mmu.142	NM_001195533	HSP90AB1	Heat shock 90kDa protein 1, beta
C06	Mmu.12316	XM_001097691	LOC709186	Hypoxanthine-guanine phosphoribosyltransferase-like
C07	Mmu.4112	XM_001102090	LOC712934	Hypothetical LOC712934
C08	Mmu.2957	XM_001091198	NONO	Non-POU domain containing, octamer-binding
C09	Mmu.4126	NM_001195658	PGK1	Phosphoglycerate kinase 1
C10	Mmu.17244	XM_001087648	PPIH	Peptidylprolyl isomerase H (cyclophilin H)
C11	Mmu.2512	XM_001115079	RPL13A	Ribosomal protein L13A
C12	Mmu.12238	XM_001085368	TBP	TATA box binding protein
D01	Mmu.4974	NM_001033084	ACTB	Actin, beta
D02	Mmu.5037	NM_001047137	B2M	Beta-2-microglobulin
D03	Mmu.3145	XM_001105471	GAPDH	Glyceraldehyde-3-phosphate dehydrogenase
D04	Mmu.9544	XM_001087699	GUSBB	Glucuronidase, beta, b
D05	Mmu.142	NM_001195533	HSP90AB1	Heat shock 90kDa protein 1, beta
D06	Mmu.12316	XM_001097691	LOC709186	Hypoxanthine-guanine phosphoribosyltransferase-like
D07	Mmu.4112	XM_001102090	LOC712934	Hypothetical LOC712934
D08	Mmu.2957	XM_001091198	NONO	Non-POU domain containing, octamer-binding
D09	Mmu.4126	NM_001195658	PGK1	Phosphoglycerate kinase 1
D10	Mmu.17244	XM_001087648	PPIH	Peptidylprolyl isomerase H (cyclophilin H)
D11	Mmu.2512	XM_001115079	RPL13A	Ribosomal protein L13A
D12	Mmu.12238	XM_001085368	TBP	TATA box binding protein
E01	Mmu.4974	NM_001033084	ACTB	Actin, beta
E02	Mmu.5037	NM_001047137	B2M	Beta-2-microglobulin
E03	Mmu.3145	XM_001105471	GAPDH	Glyceraldehyde-3-phosphate dehydrogenase
E04	Mmu.9544	XM_001087699	GUSBB	Glucuronidase, beta, b
E05	Mmu.142	NM_001195533	HSP90AB1	Heat shock 90kDa protein 1, beta
E06	Mmu.12316	XM_001097691	LOC709186	Hypoxanthine-guanine phosphoribosyltransferase-like
E07	Mmu.4112	XM_001102090	LOC712934	Hypothetical LOC712934
E08	Mmu.2957	XM_001091198	NONO	Non-POU domain containing, octamer-binding
E09	Mmu.4126	NM_001195658	PGK1	Phosphoglycerate kinase 1

Position	UniGene	GenBank	Symbol	Description
E10	Mmu.17244	XM_001087648	PPIH	Peptidylprolyl isomerase H (cyclophilin H)
E11	Mmu.2512	XM_001115079	RPL13A	Ribosomal protein L13A
E12	Mmu.12238	XM_001085368	TBP	TATA box binding protein
F01	Mmu.4974	NM_001033084	ACTB	Actin, beta
F02	Mmu.5037	NM_001047137	B2M	Beta-2-microglobulin
F03	Mmu.3145	XM_001105471	GAPDH	Glyceraldehyde-3-phosphate dehydrogenase
F04	Mmu.9544	XM_001087699	GUSBB	Glucuronidase, beta, b
F05	Mmu.142	NM_001195533	HSP90AB1	Heat shock 90kDa protein 1, beta
F06	Mmu.12316	XM_001097691	LOC709186	Hypoxanthine-guanine phosphoribosyltransferase-like
F07	Mmu.4112	XM_001102090	LOC712934	Hypothetical LOC712934
F08	Mmu.2957	XM_001091198	NONO	Non-POU domain containing, octamer-binding
F09	Mmu.4126	NM_001195658	PGK1	Phosphoglycerate kinase 1
F10	Mmu.17244	XM_001087648	PPIH	Peptidylprolyl isomerase H (cyclophilin H)
F11	Mmu.2512	XM_001115079	RPL13A	Ribosomal protein L13A
F12	Mmu.12238	XM_001085368	TBP	TATA box binding protein
G01	Mmu.4974	NM_001033084	ACTB	Actin, beta
G02	Mmu.5037	NM_001047137	B2M	Beta-2-microglobulin
G03	Mmu.3145	XM_001105471	GAPDH	Glyceraldehyde-3-phosphate dehydrogenase
G04	Mmu.9544	XM_001087699	GUSBB	Glucuronidase, beta, b
G05	Mmu.142	NM_001195533	HSP90AB1	Heat shock 90kDa protein 1, beta
G06	Mmu.12316	XM_001097691	LOC709186	Hypoxanthine-guanine phosphoribosyltransferase-like
G07	Mmu.4112	XM_001102090	LOC712934	Hypothetical LOC712934
G08	Mmu.2957	XM_001091198	NONO	Non-POU domain containing, octamer-binding
G09	Mmu.4126	NM_001195658	PGK1	Phosphoglycerate kinase 1
G10	Mmu.17244	XM_001087648	PPIH	Peptidylprolyl isomerase H (cyclophilin H)
G11	Mmu.2512	XM_001115079	RPL13A	Ribosomal protein L13A
G12	Mmu.12238	XM_001085368	TBP	TATA box binding protein
H01	Mmu.4974	NM_001033084	ACTB	Actin, beta
H02	Mmu.5037	NM_001047137	B2M	Beta-2-microglobulin
H03	Mmu.3145	XM_001105471	GAPDH	Glyceraldehyde-3-phosphate dehydrogenase
H04	Mmu.9544	XM_001087699	GUSBB	Glucuronidase, beta, b
H05	Mmu.142	NM_001195533	HSP90AB1	Heat shock 90kDa protein 1, beta
H06	Mmu.12316	XM_001097691	LOC709186	Hypoxanthine-guanine phosphoribosyltransferase-like
H07	Mmu.4112	XM_001102090	LOC712934	Hypothetical LOC712934
H08	Mmu.2957	XM_001091198	NONO	Non-POU domain containing, octamer-binding
H09	Mmu.4126	NM_001195658	PGK1	Phosphoglycerate kinase 1
H10	Mmu.17244	XM_001087648	PPIH	Peptidylprolyl isomerase H (cyclophilin H)
H11	Mmu.2512	XM_001115079	RPL13A	Ribosomal protein L13A
H12	Mmu.12238	XM_001085368	TBP	TATA box binding protein

## Related products

For optimal performance, RT<sup>2</sup> Profiler PCR Arrays should be used together with the RT<sup>2</sup> First Strand Kit for cDNA synthesis and RT<sup>2</sup> SYBR<sup>®</sup> Green qPCR Mastermixes for PCR.

Product	Contents	Cat. no.
RT <sup>2</sup> First Strand Kit (12)	Enzymes and reagents for cDNA synthesis	330401
RT <sup>2</sup> SYBR Green ROX <sup>™</sup> FAST Mastermix (2)*	For 2 x 96 assays in 96-well plates; suitable for use with the Rotor-Gene Q and other Rotor-Gene cyclers	330620

\* Larger kit sizes available; please inquire.

RT<sup>2</sup> Profiler PCR Array products are intended for molecular biology applications. These products are not intended for the diagnosis, prevention, or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at [www.qiagen.com](http://www.qiagen.com) or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN<sup>®</sup>, Rotor-Gene<sup>®</sup>, Rotor-Disc<sup>™</sup> (QIAGEN Group); ROX<sup>™</sup> (Applied Biosystems or its subsidiaries); SYBR<sup>®</sup> (Molecular Probes, Inc.).

1067688 03/2011 © 2011 QIAGEN, all rights reserved.

[www.qiagen.com](http://www.qiagen.com)

Canada ■ 800-572-9613

China ■ 8621-3865-3865

Denmark ■ 80-885945

Finland ■ 0800-914416

France ■ 01-60-920-930

Germany ■ 02103-29-12000

Hong Kong ■ 800 933 965

Ireland ■ 1800 555 049

Italy ■ 800-787980

Japan ■ 03-6890-7300

Korea (South) ■ 080-000-7145

Luxembourg ■ 8002 2076

Mexico ■ 01-800-7742-436

The Netherlands ■ 0800 0229592

Norway ■ 800-18859

Singapore ■ 1800-742-4368

Spain ■ 91-630-7050

Sweden ■ 020-790282

Switzerland ■ 055-254-22-11

UK ■ 01293-422-911

USA ■ 800-426-8157

Australia ■ 1-800-243-800

Austria ■ 0800/281010

Belgium ■ 0800-79612

Brazil ■ 0800-557779



Sample & Assay Technologies