

QuantiNova® LNA® Probe PCR Focus Panels (96-Well Format and 384-Well [4 x 96] Format)

Human EGF / PDGF Signaling Pathway

Cat. no. 249955 UPHS-040ZA

For study focus gene expression analysis

Shipping and storage

QuantiNova LNA Probe PCR Focus Panels are shipped at room temperature. Immediately upon receipt, they should be stored protected from light at 2–8°C for short term storage or at –30°C to –15°C for long time storage. Under these conditions, all components are stable for at least 12 months.

Note: Open the package and store the products appropriately immediately upon receipt.

For optimal performance, QuantiNova LNA Probe PCR Focus Panels should be used together with the QuantiNova Reverse Transcription Kit for cDNA synthesis and the QuantiNova Probe PCR Kit (Mastermix) for PCR.

Panel layout (96-well): QuantiNova LNA Probe PCR Focus Panel

For the 384-well (4 × 96) PCR panels, genes are present in a staggered format. Refer to the QuantiNova LNA Probe PCR Handbook at www.qiagen.com for further details.

	1	2	3	4	5	6	7	8	9	10	11	12
A	ACTR2	AKT1	AKT2	AKT3	ARAF	ATF1	ATF2	BAD	BCAR1	BCL2	BRAF	CASP3
B	CASP9	CBL	CCND1	CHUK	COL1A1	CREB1	CSNK2A1	CSNK2B	DUSP1	DUSP6	EGF	EGFR
C	EGR1	EIF4E	ELK1	EPS8	FASLG	FN1	FOS	FOXO3	GAB1	GRB2	GSK3A	GSK3B
D	HBEGF	HRAS	IKKKB	IL2	JAK1	JUN	KRAS	LTA	MAP2K1	MAP2K4	MAP2K7	MAP3K2
E	MAPK1	MAPK10	MAPK3	MAPK8	MAPK9	MKNK1	MMP7	NCK2	NFATC3	NFKB1	NRAS	NUF62
F	PDGFA	PDGFB	PDGFRA	PDPK1	PIK3CA	PIK3R1	PIK3R2	PLAT	PLCG1	PPP2CA	PRKCA	PTEN
G	RAF1	RAP1A	RASA1	RHOA	RPS6KA5	RPS6KB1	SHC1	SRC	STAT1	STAT3	STAT5A	TP53
H	ACTB	B2M	GAPDH	HPRT1	RPLP0	HGDC	QIC	QIC	QIC	PPC	PPC	PPC

Gene table: QuantiNova LNA Probe PCR Focus Panel

Position	Assay	Name	Symbol	Ensembl ID	Description
A01	UPFH0558245	ENST00000377982.8	ACTR2	ENSG00000138071	ARP2 actin related protein 2 homolog Source HGNC Symbol Acc HGNC 169
A02	UPFH0453992	ENST00000555528.5	AKT1	ENSG00000142208	AKT serine/threonine kinase 1 Source HGNC Symbol Acc HGNC 391
A03	UPFH0246744	ENST00000441941.6	AKT2	ENSG00000105221	AKT serine/threonine kinase 2 Source HGNC Symbol Acc HGNC 392
A04	UPFH0198803	ENST00000336199.9	AKT3	ENSG00000117020	AKT serine/threonine kinase 3 Source HGNC Symbol Acc HGNC 393
A05	UPFH0579274	ENST00000290277.10	ARAF	ENSG00000078061	A-Raf proto-oncogene, serine/threonine kinase Source HGNC Symbol Acc HGNC 646
A06	UPFH0059653	ENST00000551831.5	ATF1	ENSG00000123268	activating transcription factor 1 Source HGNC Symbol Acc HGNC 783
A07	UPFH1132243	ENST00000426833.7	ATF2	ENSG00000115966	activating transcription factor 2 Source HGNC Symbol Acc HGNC 784
A08	UPFH0437748	ENST00000394532.7	BAD	ENSG00000002330	BCL2 associated agonist of cell death Source HGNC Symbol Acc HGNC 936
A09	UPFH0300381	ENST00000418647.7	BCAR1	ENSG00000050820	BCAR1, Cas family scaffold protein Source HGNC Symbol Acc HGNC 971
A10	UPFH1132900	ENST00000333681.5	BCL2	ENSG000000171791	BCL2, apoptosis regulator Source HGNC Symbol Acc HGNC 990
A11	UPFH0559252	ENST00000646891.1	BRAF	ENSG00000157764	B-Raf proto-oncogene, serine/threonine kinase Source HGNC Symbol Acc HGNC 1097
A12	UPFH1132892	ENST00000523916.5	CASP3	ENSG00000164305	caspase 3 Source HGNC Symbol Acc HGNC 1504
B01	UPFH0083593	ENST00000440484.1	CASP9	ENSG00000132906	caspase 9 Source HGNC Symbol Acc HGNC 1511
B02	UPFH0561363	ENST00000634840.1	CBL	ENSG00000110395	Cbl proto-oncogene Source HGNC Symbol Acc HGNC 1541
B03	UPFH0430337	ENST00000227507.2	CCND1	ENSG00000110092	cyclin D1 Source HGNC Symbol Acc HGNC 1582
B04	UPFH1132316	ENST00000370397.8	CHUK	ENSG000000213341	conserved helix-loop-helix ubiquitous kinase Source HGNC Symbol Acc HGNC 1974
B05	UPFH0361104	ENST00000225964.9	COL1A1	ENSG00000108821	collagen type I alpha 1 chain Source HGNC Symbol Acc HGNC 2197
B06	UPFH0199960	ENST00000480189.5	CREB1	ENSG00000118260	cAMP responsive element binding protein 1 Source HGNC Symbol Acc HGNC 2345
B07	UPFH0055149	ENST00000645249.1	CSNK2A1	ENSG00000101266	casein kinase 2 alpha 1 Source HGNC Symbol Acc HGNC 2457
B08	UPFH0261524	ENST00000375882.6	CSNK2B	ENSG000000204435	casein kinase 2 beta Source HGNC Symbol Acc HGNC 2460
B09	UPFH1132372	ENST00000239223.4	DUSP1	ENSG00000120129	dual specificity phosphatase 1 Source HGNC Symbol Acc HGNC 3064
B10	UPFH0610812	ENST00000547291.1	DUSP6	ENSG00000139318	dual specificity phosphatase 6 Source HGNC Symbol Acc HGNC 3072
B11	UPFH1132380	ENST00000503392.1	EGF	ENSG00000138798	epidermal growth factor Source HGNC Symbol Acc HGNC 3229
B12	UPFH1132381	ENST00000420316.6	EGFR	ENSG00000146648	epidermal growth factor receptor Source HGNC Symbol Acc HGNC 3236
C01	UPFH0558832	ENST00000239938.5	EGR1	ENSG00000120738	early growth response 1 Source HGNC Symbol Acc HGNC 3238
C02	UPFH0008084	ENST00000515638.5	EIF4E	ENSG00000151247	eukaryotic translation initiation factor 4E Source HGNC Symbol Acc HGNC 3287
C03	UPFH0614738	ENST00000376983.8	ELK1	ENSG00000126767	ELK1, ETS transcription factor Source HGNC Symbol Acc HGNC 3321
C04	UPFH0507645	ENST00000642939.1	EPS8	ENSG00000151491	epidermal growth factor receptor pathway substrate 8 Source HGNC Symbol Acc HGNC 3420
C05	UPFH1132396	ENST00000367721.3	FASLG	ENSG00000117560	Fas ligand Source HGNC Symbol Acc HGNC 11936
C06	UPFH0605066	ENST00000336916.8	FN1	ENSG00000115414	fibronectin 1 Source HGNC Symbol Acc HGNC 3778
C07	UPFH1132401	ENST00000555242.1	FOS	ENSG00000170345	Fos proto-oncogene, AP-1 transcription factor subunit Source HGNC Symbol Acc HGNC 3796
C08	UPFH0442611	ENST00000540898.1	FOXO3	ENSG00000118689	forkhead box O3 Source HGNC Symbol Acc HGNC 3821
C09	UPFH0243246	ENST00000509992.1	GAB1	ENSG00000109458	GRB2 associated binding protein 1 Source HGNC Symbol Acc HGNC 4066
C10	UPFH1132426	ENST00000392564.5	GRB2	ENSG00000177885	growth factor receptor bound protein 2 Source HGNC Symbol Acc HGNC 4566
		ENST00000398		ENSG000000	

Position	Assay	Name	Symbol	Ensembl ID	Description
C11	UPFH1132428	249.8	GSK3A	105723	glycogen synthase kinase 3 alpha Source HGNC Symbol Acc HGNC 4616
C12	UPFH0470775	ENST00000316626.5	GSK3B	ENSG00000082701	glycogen synthase kinase 3 beta Source HGNC Symbol Acc HGNC 4617
D01	UPFH0129784	ENST00000230990.7	HBEGF	ENSG00000113070	heparin binding EGF like growth factor Source HGNC Symbol Acc HGNC 3059
D02	UPFH1132981	ENST00000311189.8	HRAS	ENSG00000174775	HRas proto-oncogene, GTPase Source HGNC Symbol Acc HGNC 5173
D03	UPFH0596293	ENST00000649612.2	IKKB	ENSG00000104365	inhibitor of nuclear factor kappa B kinase subunit beta Source HGNC Symbol Acc HGNC 5960
D04	UPFH0116492	ENST00000226730.4	IL2	ENSG00000109471	interleukin 2 Source HGNC Symbol Acc HGNC 6001
D05	UPFH1132963	ENST00000342505.5	JAK1	ENSG00000162434	Janus kinase 1 Source HGNC Symbol Acc HGNC 6190
D06	UPFH0569765	ENST00000371222.3	JUN	ENSG00000177606	Jun proto-oncogene, AP-1 transcription factor subunit Source HGNC Symbol Acc HGNC 6204
D07	UPFH0376060	ENST00000557334.5	KRAS	ENSG00000133703	KRAS proto-oncogene, GTPase Source HGNC Symbol Acc HGNC 6407
D08	UPFH1132824	ENST00000454783.5	LTA	ENSG00000226979	lymphotoxin alpha Source HGNC Symbol Acc HGNC 6709
D09	UPFH1132932	ENST00000307102.9	MAP2K1	ENSG00000169032	mitogen-activated protein kinase kinase 1 Source HGNC Symbol Acc HGNC 6840
D10	UPFH1132531	ENST00000353533.10	MAP2K4	ENSG00000065559	mitogen-activated protein kinase kinase 4 Source HGNC Symbol Acc HGNC 6844
D11	UPFH1124074	ENST00000468058.1	MAP2K7	ENSG00000076984	mitogen-activated protein kinase kinase 7 Source HGNC Symbol Acc HGNC 6847
D12	UPFH0130844	ENST00000409947.5	MAP3K2	ENSG00000169967	mitogen-activated protein kinase kinase kinase 2 Source HGNC Symbol Acc HGNC 6854
E01	UPFH0366815	ENST00000215832.10	MAPK1	ENSG00000100030	mitogen-activated protein kinase 1 Source HGNC Symbol Acc HGNC 6871
E02	UPFH0530550	ENST00000641563.1	MAPK10	ENSG00000109339	mitogen-activated protein kinase 10 Source HGNC Symbol Acc HGNC 6872
E03	UPFH1132534	ENST00000481230.1	MAPK3	ENSG00000102882	mitogen-activated protein kinase 3 Source HGNC Symbol Acc HGNC 6877
E04	UPFH1132535	ENST00000374179.8	MAPK8	ENSG00000107643	mitogen-activated protein kinase 8 Source HGNC Symbol Acc HGNC 6881
E05	UPFH0148170	ENST00000397072.7	MAPK9	ENSG00000050748	mitogen-activated protein kinase 9 Source HGNC Symbol Acc HGNC 6886
E06	UPFH0234693	ENST0000049619.6	MKNIK1	ENSG00000079277	MAP kinase interacting serine/threonine kinase 1 Source HGNC Symbol Acc HGNC 7110
E07	UPFH0230006	ENST00000260227.5	MMP7	ENSG00000137673	matrix metalloproteinase 7 Source HGNC Symbol Acc HGNC 7174
E08	UPFH0049450	ENST00000425756.1	NCK2	ENSG00000071051	NCK adaptor protein 2 Source HGNC Symbol Acc HGNC 7665
E09	UPFH0123113	ENST00000562926.5	NFATC3	ENSG00000072736	nuclear factor of activated T cells 3 Source HGNC Symbol Acc HGNC 7777
E10	UPFH1132828	ENST00000226574.9	NFKB1	ENSG00000109320	nuclear factor kappa B subunit 1 Source HGNC Symbol Acc HGNC 7794
E11	UPFH0189023	ENST00000369535.5	NRAS	ENSG00000213281	NRAS proto-oncogene, GTPase Source HGNC Symbol Acc HGNC 7989
E12	UPFH0319644	ENST00000596217.1	NUP62	ENSG00000213024	nucleoporin 62 Source HGNC Symbol Acc HGNC 8066
F01	UPFH1132608	ENST00000354513.9	PDGFA	ENSG00000197461	platelet derived growth factor subunit A Source HGNC Symbol Acc HGNC 8799
F02	UPFH0524988	ENST00000381551.8	PDGFB	ENSG00000100311	platelet derived growth factor subunit B Source HGNC Symbol Acc HGNC 8800
F03	UPFH1132609	ENST00000508170.5	PDGFRA	ENSG00000134853	platelet derived growth factor receptor alpha Source HGNC Symbol Acc HGNC 8803
F04	UPFH0104991	ENST00000389224.7	PDPK1	ENSG00000140992	3-phosphoinositide dependent protein kinase 1 Source HGNC Symbol Acc HGNC 8816
F05	UPFH0109251	ENST00000462255.1	PIK3CA	ENSG00000121879	phosphatidylinositol-4,5-bisphosphate 3-kinase catalytic subunit alpha Source HGNC Symbol Acc HGNC 8975
F06	UPFH0549122	ENST00000521381.6	PIK3R1	ENSG00000145675	phosphoinositide-3-kinase regulatory subunit 1 Source HGNC Symbol Acc HGNC 8979
F07	UPFH0463753	ENST00000617130.4	PIK3R2	ENSG00000105647	phosphoinositide-3-kinase regulatory subunit 2 Source HGNC Symbol Acc HGNC 8980
F08	UPFH1132830	ENST00000352041.7	PLAT	ENSG00000104368	plasminogen activator, tissue type Source HGNC Symbol Acc HGNC 9051
F09	UPFH0088354	ENST00000617873.4	PLCG1	ENSG00000124181	phospholipase C gamma 1 Source HGNC Symbol Acc HGNC 9065
F10	UPFH0120233	ENST00000522385.1	PPP2CA	ENSG00000113575	protein phosphatase 2 catalytic subunit alpha Source HGNC Symbol Acc HGNC 9299

Position	Assay	Name	Symbol	Ensembl ID	Description
F11	UPFH0607768	ENST00000578063.5	PRKCA	ENSG00000154229	protein kinase C alpha Source HGNC Symbol Acc HGNC 9393
F12	UPFH1132982	ENST00000371953.8	PTEN	ENSG00000171862	phosphatase and tensin homolog Source HGNC Symbol Acc HGNC 9588
G01	UPFH0380839	ENST00000416093.1	RAF1	ENSG00000132155	Raf-1 proto-oncogene, serine/threonine kinase Source HGNC Symbol Acc HGNC 9829
G02	UPFH0233643	ENST00000369709.3	RAP1A	ENSG00000116473	RAP1A, member of RAS oncogene family Source HGNC Symbol Acc HGNC 9855
G03	UPFH0086954	ENST00000509953.1	RASA1	ENSG00000145715	RAS p21 protein activator 1 Source HGNC Symbol Acc HGNC 9871
G04	UPFH1132657	ENST00000445425.4	RHOA	ENSG00000067560	ras homolog family member A Source HGNC Symbol Acc HGNC 667
G05	UPFH1122533	ENST00000556178.5	RPS6KA5	ENSG00000100784	ribosomal protein S6 kinase A5 Source HGNC Symbol Acc HGNC 10434
G06	UPFH1132665	ENST00000406116.7	RPS6KB1	ENSG00000108443	ribosomal protein S6 kinase B1 Source HGNC Symbol Acc HGNC 10436
G07	UPFH0141239	ENST00000368449.8	SHC1	ENSG00000160691	SHC adaptor protein 1 Source HGNC Symbol Acc HGNC 10840
G08	UPFH0308412	ENST00000472968.1	SRC	ENSG00000197122	SRC proto-oncogene, non-receptor tyrosine kinase Source HGNC Symbol Acc HGNC 11283
G09	UPFH1132696	ENST00000392323.6	STAT1	ENSG00000115415	signal transducer and activator of transcription 1 Source HGNC Symbol Acc HGNC 11362
G10	UPFH0531262	ENST00000404395.3	STAT3	ENSG00000168610	signal transducer and activator of transcription 3 Source HGNC Symbol Acc HGNC 11364
G11	UPFH0458805	ENST00000591556.1	STAT5A	ENSG00000126561	signal transducer and activator of transcription 5A Source HGNC Symbol Acc HGNC 11366
G12	UPFH0565795	ENST00000269305.8	TP53	ENSG00000141510	tumor protein p53 Source HGNC Symbol Acc HGNC 11998
H01	UPFH1132936	ENST00000646664.1	ACTB	ENSG00000075624	actin beta Source HGNC Symbol Acc HGNC 132
H02	UPFH1132937	ENST00000544417.5	B2M	ENSG00000166710	beta-2-microglobulin Source HGNC Symbol Acc HGNC 914
H03	UPFH1132938	ENST00000229239.10	GAPDH	ENSG00000111640	glyceraldehyde-3-phosphate dehydrogenase Source HGNC Symbol Acc HGNC 4141
H04	UPFH1132939	ENST00000298556.8	HPRT1	ENSG00000165704	hypoxanthine phosphoribosyltransferase 1 Source HGNC Symbol Acc HGNC 5157
H05	UPFH1132941	ENST00000392514.9	RPLP0	ENSG00000089157	ribosomal protein lateral stalk subunit P0 Source HGNC Symbol Acc HGNC 10371
H06	UPFH1126608	UPL_HGDC	HGDC	UPL_HGDC	Human Genomic DNA Contamination
H07	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H08	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H09	UPFH1126606	UPL_QIC	QIC	UPL_QIC	QuantiNova Internal Control
H10	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control
H11	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control
H12	UPFH1126605	UPL_PPC	PPC	UPL_PPC	Positive PCR Control



Related products

Product	Contents	Cat. no.
QuantiNova LNA Probe PCR QC Panel	These panels are designed to assess the quality of RNA samples before characterization using QuantiNova LNA Probe PCR Focus Panels; available in 96-well, 384-well, and Rotor-Disc 100 formats	249945
QuantiNova Reverse Transcription Kit (10)*	For 10 x 20 μ l reactions: 20 μ l 8x gDNA Removal Mix, 10 μ l Reverse Transcription Enzyme, 40 μ l Reverse Transcription Mix (containing RT primers), 20 μ l Internal Control RNA, 1.9 ml RNase-Free Water	205410
QuantiNova Probe RT-PCR Kit (100)*	For 100 x 20 μ l reactions: 1 ml QuantiNova Probe RT-PCR Master Mix, 20 μ l QuantiNova Probe RT Mix, 20 μ l Internal Control RNA, 500 μ l Yellow Template Dilution Buffer, 250 μ l ROX Reference Dye, 1.9 μ l RNase-Free Water	208352
QuantiNova Probe PCR Kit (100)*	For 100 x 20 μ l reactions: 1 ml 2x QuantiNova Probe PCR Master Mix, 500 μ l QuantiNova Yellow Template Dilution Buffer, 250 μ l QN ROX Reference Dye, 1.9 ml Water	208252

*Larger kit sizes available.

The QuantiNova LNA Probe PCR Focus Panels are intended for molecular biology applications. These products are not intended for the diagnosis, prevention or treatment of a disease.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, LNA®, QuantiNova®, Sample to Insight® (QIAGEN Group); SYBR® (Life Technologies Corp.). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

09/2019 © 2019 QIAGEN, all rights reserved.