

miRCURY LNA miRNA miRNome PCR Panels

Rhesus Macaque Panel 1+2 Product Data Sheet

Cat. no. 339322 YAQQ-062Z

For mature miRNA expression profiling using real-time PCR

| Format | Suitable real-time cyclers | Plate | Cat. no. |
|--------|---|---------|-------------|
| A | Applied Biosystems® models 5700, 7000, 7300, 7500, 7700, 7900HT, ViiA™ 7 (96-well block); Bio-Rad® models iCycler®, iQ™5, MyiQ™, MyiQ2; Bio-Rad/MJ Research Chromo4™; Eppendorf® MasterCycler® ep realplex models 2, 2s, 4, 4s; Stratagene® models Mx3005P®, Mx3000P®; Takara: TP-800 | 96-well | YAQQ-062Z A |
| C | Applied Biosystems models 7500 (Fast block), 7900HT (Fast block), StepOnePlus™, ViiA 7 (Fast block) | 96-well | YAQQ-062Z C |
| D | Bio-Rad CFX96™; Bio-Rad/MJ Research models DNA Engine Opticon®, DNA Engine Opticon 2; Stratagene Mx4000® | 96-well | YAQQ-062Z D |
| F | Roche® LightCycler® 480 (96-well block) | 96-well | YAQQ-062Z F |

Description

The miRCURY LNA miRNA miRNome PCR Panels profile the expression of hundreds of the most abundantly expressed and best characterized miRNA sequences in the a given species-specific miRNome as annotated in miRBase (www.mirbase.org). Although they are well characterized, each of these miRNA sequences can regulate one to several messenger RNA transcripts, and conversely one mRNA can be regulated by one to several miRNA sequences. Therefore, the complex role that any given known miRNA sequence plays has yet to be completely defined. Use of the miRCURY LNA miRNA miRNome PCR Panels maximizes the likelihood of discovering miRNA sequences whose expression patterns correlate with the biological phenotypes under study. A set of controls included on each plate enables data analysis using the $\Delta\Delta CT$ method of relative quantification, assessment of reverse transcription performance, and assessment of PCR performance. The Human miRNome miScript miRNA PCR Array enables easy and reliable, SYBR Green-based real-time PCR analysis of the expression of 1008 miRNA sequences and discovery of those most important to the research area of interest.

For further details, consult the *miRCURY LNA™ miRCURY LNA™ miRNA Focus PCR Panels Handbook*.

Ordering Information

| Product | Contents | Cat. no. |
|---|--|----------|
| miRCURY LNA miRNA Focus PCR Panels | miRCURY LNA miRNA PCR Panels for application-based miRNome profiling, available in 96-well or 384-well format; for SYBR® Green-based detection | 339325 |
| miRCURY LNA miRNA miRNome PCR Panels | miRCURY LNA miRNA PCR Panels for PCR-based miRNome profiling, available in 384-well format; for SYBR® Green-based detection | 339322 |
| miRCURY LNA miRNA QC PCR Panel | miRCURY LNA miRNA PCR Panel of quality control assays, available in 96-well or 384-well format; for SYBR® Green-based detection | 339331 |
| miRCURY LNA miRNA Custom PCR Panels | 8 identical, ready-to-use 96- or 384-well plates; each well contains primers sufficient for one 10 µl reaction; for SYBR® Green-based detection | 339330 |
| miRCURY LNA Custom PCR Panel Additional Plate | Additional miRCURY LNA Custom PCR Panel plates; set of 4 plates; only available in addition to the base plates ordered through the core product (cat. no. 339330) | 339332 |
| miRCURY LNA miRNA PCR Assays | Contains forward and reverse primers for 200 SYBR® Green-based, real-time qPCR reactions, 166 EvaGreen-based digital PCR reactions for Nanoplate 8.5k or 50 EvaGreen-based digital PCR reactions for Nanoplate 26k | 339306 |
| miRCURY LNA miRNA Probe PCR Assays | Complete premixed assays containing LNA-enhanced target-specific forward primer and probe. For 200 reactions. | 339350 |
| miRCURY LNA miRNA Custom Probe PCR Assays | Custom-designed, target-specific forward primer and probe for any user-defined miRNA target. Complete premixed assay for 200 reactions. | 339351 |

Related Products

| Product | Contents | Cat. no. |
|---|--|----------|
| miRCURY LNA RT Kit | For 8–64 cDNA synthesis reactions: 5x RT SYBR Green Reaction Buffer, 5x RT Probe Reaction Buffer, 10x RT Enzyme Mix, UniSp6, RNA Spike-in template, RNase-Free Water | 339340 |
| RNA Spike-In Kit, For RT | Contains the UniSp2, UniSp4, and UniSp5 RNA Spike-in Template Mix and the cel-miR-39-3p RNA Spike-in Template | 339390 |
| miRCURY LNA SYBR® Green PCR Kits (200) | For 200 reactions: 2X miRCURY SYBR Green Master Mix, RNase-Free Water | 339345 |
| miRCURY LNA SYBR® Green PCR Kits (600) | For 600 reactions: 2X miRCURY SYBR Green Master Mix, RNase-Free Water | 339346 |
| miRCURY LNA SYBR® Green PCR Kits (4000) | For 4000 reactions: 2X miRCURY SYBR Green Master Mix, RNase-Free Water | 339347 |
| miRCURY LNA Probe PCR Kit (200) | For 200 reactions: 2X QuantiNova Probe Master Mix, 10X miRCURY Probe Univ. Primer, Rox Reference Dye, RNase-Free Water | 339371 |
| miRCURY LNA Probe PCR Kit (800) | For 800 reactions: 2X QuantiNova Probe Master Mix, 10X miRCURY Probe Univ. Primer, Rox Reference Dye, RNase-Free Water | 339372 |
| miRCURY LNA Probe PCR Kit (4000) | For 4000 reactions: 2X QuantiNova Probe Master Mix, 10X miRCURY Probe Univ. Primer, Rox Reference Dye, RNase-Free Water | 339373 |
| miRCURY LNA miRNA PCR Starter Kit | Two miRCURY LNA PCR Assays of your choice, spike-in control Assay (UniSp6), one candidate endogenous control assay (miR-103-3p) and all reagents for 20 reverse transcription reactions and 100 PCR amplifications; for SYBR® Green-based qPCR detection | 339320 |

Important

The miRCURY LNA™ miRNA PCR Panels are Ready-to-Use and designed for optimal performance with the miRCURY LNA RT Kit and the miRCURY LNA SYBR® Green PCR Kit. The performance of the primer sets will be affected when used in combination with less than optimal reagents. Use the miRCURY LNA miRNA PCR Panel Handbook for experiment setup. RNA work requires specific handling and precautions should be taken to prevent RNase contamination and degradation of the RNA sample and reagents.

Shipping and storage

The Ready-to-Use miRCURY LNA miRNA PCR Panels are shipped at room temperature and can be stored at 4°C for at least 6 months. For long term storage, it is recommended to place the panels at -20°C. Under these conditions, the LNA PCR primers are stable for at least 6 months after receipt.

Intended use

miRCURY LNA miRNA PCR Assays are intended for molecular biology applications. This product is not intended for the diagnosis, prevention or treatment of a disease.

Safety information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. For more information, please consult the appropriate safety data sheets (SDSs). These are available online in convenient and compact PDF format at **www.qiagen.com/safety** where you can find, view and print the SDS for each QIAGEN kit and kit component.

For up-to-date licensing information and product-specific disclaimers, see the respective QIAGEN kit handbook or user manual. QIAGEN kit handbooks and user manuals are available at www.qiagen.com or can be requested from QIAGEN Technical Services or your local distributor.

Trademarks: QIAGEN®, Sample to Insight®, GeneGlobe®, LNA™ (QIAGEN Group), SYBR® (Life Technologies Corporation). Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law. 1107212 05/2017 HB-2376-001 © 2017 QIAGEN, all rights reserved.
